

D20 II

The compact basic solution for the professional use in small labs and dental clinics.

Rapid Shape D20 II produces high-precision and esthetic dental products made from certified, biocompatible materials.

Thanks to short production times, there is no need for patients to schedule a second appointment. The D20 II is the perfect compact unit for medical practices, clinics, and small labs. A certified software interface for common lab programs comes standard in each unit.

Performance parameters	D20 II
Building area	130 × 75 mm
Native pixel	+/- 34 µm
Max. part height	90 mm
Light source	385 nm UV LED
Resolution	HD 1920 × 1080 px
Dimension (W × H × D)	335 × 541 × 349 mm
Connections	WLAN, TCP/IP, USB
Control	7" LCD-Display, 3-button control

